


ENTERTAINMENT SYSTEM MOUNTING KIT MEDICAL GRADE & STANDARD TV SYSTEMS

*TV not included with kit;
for display purposes only



Many of our hospital and clinic partners would like to utilize their own TV with the Sechrist Entertainment System. This option is available and made convenient through the Sechrist Entertainment Mounting Kit which includes all the necessary parts and equipment to accommodate both Medical Grade TV Systems and Standard TV Systems (power cord sold separately). Easy to follow instructions allow for simple installation by either a Sechrist Authorized Technician or hospital/clinic biomed personnel. The Mounting Kit provides flexibility when selecting the Entertainment System for your hyperbaric patients.



*TV not included with kit; for display purposes only

THE WORLD LEADER IN HYPERBARIC TECHNOLOGY & MEDICINE

ENTERTAINMENT SYSTEM MOUNTING KIT

MEDICAL GRADE & STANDARD TV SYSTEMS

MOUNTING KITS	
FOR STANDARD TV SYSTEM	PART NUMBER
Models 3200 / 3300E / 3300H / 3600E / 3600H Series Hyperbaric Chambers	21868
Model 4100H Series Hyperbaric Chambers	21866
FOR MEDICAL GRADE TV SYSTEM	
Models 3200 / 3300E / 3300H / 3600E / 3600H Series Hyperbaric Chambers	21869
Model 4100H Series Hyperbaric Chambers	21867
TECHNICAL RECOMMENDATIONS PER SYSTEM	
TV SIZE	TV must not weigh more than 26 pounds
MOUNTING PATTERN	TV mounting pattern must be VESA 75x75mm or 100x100mm. The sum of TV thickness and distance from TV mount center to TV bottom must not exceed 13 inches.
AUDIO OUTPUT STANDARD TV	Audio output must have 3.5mm or 1/4 in. headphone jack output if utilizing Sechrist audio output cable
AUDIO OUTPUT MEDICAL GRADE TV	Audio output must have 6-pin pillow port output if utilizing Sechrist audio output cable
POWER CORD	TV power cord must be sufficiently long to reach from outlet to top of chamber**, 12 ft. recommended
KIT COMPONENTS AND PART NUMBERS	
AUDIO OUTPUT CABLE FOR STANDARD TV	TV audio output standard cable 3.5mm headphone jack to male dual RCA; (P/N 22277) Qty 1, 15ft
AUDIO OUTPUT CABLE FOR MEDICAL GRADE TV	TV audio output cable 6-pin pillow port to male single RCA; (P/N 22248) Qty 1, 12 ft.
AUDIO CABLE	Audio cable: 1/4 in. phone plug to male dual RCA; (P/N 22247) Qty 1, 12 ft.
CABLES & HARDWARE	Coaxial Cable, Video low profile; (P/N 22245) Qty 1, 25 ft. Tubing Split Corrugated 3/4 in. - inner diameter (P/N 43541) Qty 1, 12 ft. Cable Holder Adhesive back (P/N 43715) Qty 5 Ground Loop Isolator (P/N 41063) Qty 1 Screw, Panhead Phillips M4 x 6 mm Long (P/N 56541) and 8mm Long (P/N 46541) Qty 4 each
POWER CORD (not included in kit) **	Power Cord, 13A, 125VAC, IEC320 female receptacle; (P/N 22223), 15 ft **Note - please verify that cord receptacle fits TV purchased (see image below).



**Shown: Power Cord Receptacle
(Part number 22223)

sechristusa.com

THE WORLD LEADER IN HYPERBARIC TECHNOLOGY & MEDICINE